

Hard-to-heal venous leg ulcers: Results from a double-blind RCT of a silver foam dressing*

Patricia Senet, Hôpital Tenon, Paris, France · Birte Petersen, Coloplast, Humlebæk, Denmark

Introduction

Hard-to-heal¹ ulcers are frequently inflammatory and often highly colonized. Ionised silver has anti-inflammatory and antibacterial properties. Previous studies have shown that a silver foam dressing* provides a faster reduction in wound area of infected wounds than a dressings without silver^{1,2}.

Aim

To investigate the efficacy of a silver foam dressing* in comparison with a foam dressing without silver⁺ for the treatment of hard-to-heal venous leg ulcers.

Methods

This is a subgroup analysis of the French patient group (n=75) from a 10-week, multinational, double-blind randomised, controlled trial (RCT) of patients with a venous or predominantly venous leg ulcer that had failed to heal despite appropriate therapy, including compression therapy, in the 4 weeks prior to inclusion.

The patients were centrally randomized and allocated to treatment with the silver foam or the comparator foam for 6 weeks by an interactive voice response service (IVRS). Thereafter both groups continued with the comparator foam for 4 weeks.


The primary endpoint was relative reduction in ulcer area at week 6. Data from week 10 are also presented. Wound size and area were measured by planimetry. Data were analysed on the intention-to-treat (ITT) population with last observation carried forward (LOCF) using ANCOVA with baseline area, age and BMI as covariates and gender, ulcer age and treatment as fixed effects.

Results

Baseline parameters were comparable in the two treatment groups.

| Baseline parameters | Biatain | Biatain Ag | p-value |
|---|-------------|-------------|---------|
| N | 38 | 37 | |
| Gender | 23 /15 | 22/15 | 0.92 |
| Female / Male | | | |
| Age (years) Mean (SD) | 76.2 (9.5) | 80.4 (7.8) | 0.04 |
| Ulcer area (cm ² , planimetry) Mean (SD) | 18.5 (15.8) | 17.1 (13.7) | 0.69 |
| Ulcer duration (years) Mean (SD) | 3.1 (3.9) | 3.0 (5.3) | 0.90 |
| Ankle/Brachial/Index (ABI) Mean (SD) | 1 (0.1) | 1 (0.2) | 0.25 |
| Ulcer origin (N, %) | | | |
| Deep venous insufficiency | 24 (64.9%) | 21 (56.8%) | 0.48 |
| Superficial venous insufficiency | 16 (43.2%) | 14 (37.8%) | 0.81 |
| History of phlebitis | 14 (37.8%) | 17 (45.9%) | 0.64 |

Figure 1. After 6 weeks ulcer area was reduced by 44% with the silver foam and by 22% with the comparator (p=0.0229). After 10 weeks (when comparator foam had been applied on both groups for 4 weeks) the ulcer area was reduced by 60% in the silver foam group and by 30% in the comparator group (p=0.0262).



The linear wound healing rate (Gilman) was 0.63 mm/week for the silver foam and 0.33 mm/week for the comparator.

Two adverse events possibly related to treatment were registered. In the silver foam group, one patient had pain in the study ulcer (severe) and in the comparator group one patient had pruritus (itching) (moderate).

Conclusions

In this subgroup analysis of the French patient group from a double-blind RCT, a silver foam dressing was effective for the treatment of hard-to-heal venous leg ulcers in comparison with a similar foam dressing without silver. This effect was still significant after 10 weeks, 4 weeks after treatment was switched from the silver foam to the comparator foam.

References

1. Jørgensen et al. Int Wound J 2005;2(1):64-73.
2. Münster et al. J Wound Care 2006;15(5):199-206.

*Biatain Ag, *Biatain, Coloplast A/S

Biatain® Ag

– superior absorption for infected wounds

New evidence

Hard-to-heal venous leg ulcers: Results from a double-blind RCT of a silver foam dressing. Senet and Petersen. EWMA 2012

Superior
absorption
for faster
healing

Key findings

- Biatain Ag effectively started healing of hard-to-heal venous leg ulcers
 - ➔ Wound area was reduced by 44% at week 6
- The effect persisted even after Biatain Ag was switched to Biatain without silver after 6 weeks
 - ➔ Wound area was reduced by 60% at week 10
- Healing was twice as fast with Biatain Ag during the full 10 weeks study period

Description of Biatain Ag

Biatain Ag is a soft and conformable silver foam dressing with superior absorption that helps infected wounds heal faster^{1,2}. Major independent studies have proven that Biatain Ag is the only wound dressing that is effective on all bacteria commonly found in infected wounds^{3,4}.

References: 1. Jørgensen et al. International Wound Journal 2005;2(1):64-73 2. Münter et al. Journal of Wound Care. 2006;15(5):199-206 3. Ip et al. Journal of Medical Microbiology 2006;55:59-63 4. Basterzi et al. Wounds July 2010.

Biatain® Silicone Ag



| Size | Itemno. | National code |
|-----------|---------|---------------|
| 7.5x7.5 | 39636 | |
| 10x10 | 39637 | |
| 12.5x12.5 | 39638 | |

Biatain® Ag Non-Adhesive



| Size | Itemno. | National code |
|-------------|---------|---------------|
| 5x7 | 5105 | |
| 10x10 | 9622 | |
| 10x20 | 9623 | |
| 15x15 | 9625 | |
| 20x20 | 9626 | |
| Soft Cavity | 9628 | |

Biatain® Ag Adhesive



| Size | Itemno. | National code |
|--------------|---------|---------------|
| 7.5x7.5 | 9631 | |
| 12.5x12.5 | 9631 | |
| 15x15 | 3484 | |
| 18x18 | 9635 | |
| 17x17 Sacral | 9641 | |
| 19x20 Heel | 9643 | |

Ostomy Care
Urology & Continence Care
Wound & Skin Care



Coloplast develops products and services that make life easier for people with very personal and private medical conditions. Working closely with the people who use our products, we create solutions that are sensitive to their special needs. We call this intimate healthcare. Our business includes ostomy care, urology and continence care and wound and skin care. We operate globally and employ more than 7,000 people.

The Coloplast logo is a registered trademark of Coloplast A/S. © 2012-04.
All rights reserved Coloplast A/S, 3050 Humlebæk, Denmark.

Coloplast A/S
Holtedam 1
3050 Humlebæk
Denmark

www.coloplast.com