

Biatain[®] Alginate

- Biatain Alginate formerly known as SeaSorb[®] Soft
- A unique dressing which combines the best from alginates and hydrofibers in one
- Easy to remove in one piece




Biatain[®] Alginate

- Superior absorption for sloughy wounds and cavity filling

Biatain Alginate is a highly absorbent alginate dressing for moderately to heavily exuding wounds of any size, shape and depth. Biatain Alginate has superior absorption compared to other alginate and hydrofiber dressings².

- Ideal for cavity filling³
- Unique dressing integrity^{1,2}
- Minimal risk of leakage and maceration
- Documented haemostatic effect⁴



On what types of wounds would I use Biatain Alginate?


Biatain Alginate can be used on all moderate to heavily exuding partial to full thickness wounds. Biatain Alginate may assist in supporting the control of minor bleeding in superficial wounds. Biatain Alginate can be used with compression therapy. Biatain Alginate should not be used on dry or lightly exuding wounds or to control heavy bleeding.

Biatain Alginate is indicated for:

- Diabetic foot ulcers
- Donor sites
- Second degree burns
- Recommended wear time up to 7 days dependant on exudate levels
- Pressure ulcers
- Leg ulcers
- Sloughy wounds



Biatain Alginate Wound Dressing

	Size (cm)	Qty	Code	NHS	PIP
	5x5	30	3705	ELS068	377-6747
	10x10	10	3710	ELS069	377-6754
	15x15	10	3715	ELS070	377-6762
	3x44 filler	6	3740	ELS071	377-6770

References

1. Thomas et al. www.dressings.org/TechnicalPublications/PDF/Coloplast-Dressings-Testing-2003-2004.pdf
2. Timmons. Professional Nurse 1999;14(7).
3. <http://www.worldwidewounds.com/1998/june/Alginates-FAQ/alginate-questions.html>
4. Segal et al. Journal of Biomaterials Applications 1998;249.

Ostomy Care
Urology & Continence Care
Wound & Skin Care

Coloplast develops products and services that make life easier for people with very personal and private medical conditions. Working closely with the people who use our products, we create solutions that are sensitive to their special needs. We call this intimate healthcare. Our business includes ostomy care, urology and continence care and wound and skin care. We operate globally and employ more than 7,000 people.

The Coloplast logo is a registered trademark of Coloplast A/S. © 09.2011
All rights reserved Coloplast A/S, 3050 Humlebæk, Denmark.



Coloplast Ltd
Peterborough Business Park
Lynchwood
Peterborough
PE2 6FX

CM621N

www.coloplast.co.uk/woundandskincare