

Biatain® Super

- Biatain® Super - formerly known as Alione®
- Super absorbent dressing for use on highly exuding wounds

Biatain® Super - superior absorption for highly exuding wounds

What makes Biatain® Super highly absorbent?

Biatain Super's hydrocapillary pad contains super absorbent particles, which rapidly absorb and lock exudate away from the wound bed and surrounding skin. The semi-permeable top film allows evaporation, which increases the overall capacity of the dressing.



Benefits of Biatain Super for your patients

- Minimal maceration resulting in faster wound healing
- Reduced leakage so faster healing
- Longer wear time¹ resulting in fewer dressing changes
- Bacteria proof² and water resistant with a semi permeable layer enabling the patient to shower
- Recommended wear time up to 7 days dependent on exudate levels




On what types of wounds would I use Biatain Super?

- Biatain Super is most suitable for use on highly exuding wounds, such as: pressure ulcers, leg ulcers, non-infected diabetic foot ulcers, second degree burns and surgical wounds
- Biatain Super can be used under compression bandaging
- Biatain Super non-adhesive dressings are suitable for patients with dry or fragile skin



**Superior
Absorption™**

Biatain Super Hydrocapillary Dressing



Biatain Super Hydrocapillary Adhesive Dressing



Size (cm)	Qty	Code	NHS	PIP
10x10	10	4610	ELY103	290-2054
12.5x12.5	10	4612	ELY104	290-1999
12x20	10	4625	ELM085	302-9592
15x15	10	4615	ELY105	290-2021
20x20	10	4620	ELY114	294-1029

Biatain Super Hydrocapillary Non-Adhesive Dressing



Size (cm)	Qty	Code	NHS	PIP
10x10	10	4630	ELY106	290-1858
12.5x12.5	10	4632	ELY107	290-1957
12x20	10	4645	ELM086	302-9584
15x15	10	4635	ELY108	290-1924
20x20	10	4639	ELY115	294-1011

References

1. Haase, L. et al. Safety and performance of a new non-adhesive hydrocapillary dressing in the treatment of venous leg ulcers, presented at the 12th Conference of the European Wound Management Association.
2. S.M.T.L., Bacterial Barrier Testing Report. No: 01/1420/01. December 2001.

Ostomy Care
Urology & Continence Care
Wound & Skin Care

Coloplast develops products and services that make life easier for people with very personal and private medical conditions. Working closely with the people who use our products, we create solutions that are sensitive to their special needs. We call this intimate healthcare. Our business includes ostomy care, urology and continence care and wound and skin care. We operate globally and employ more than 7,000 people.

The Coloplast logo is a registered trademark of Coloplast A/S. © 09.2011
All rights reserved Coloplast A/S, 3050 Humlebæk, Denmark.

CM618N



Coloplast Ltd
Peterborough Business Park
Lynchwood
Peterborough
PE2 6FX

www.coloplast.co.uk/woundandskincare