

**Comfeel®**

# Ensured protection. Faster healing.

Protection and monitoring of skin and wound




Heals and Seals: Comforts and Protects

# More than 20 years of clinical experience

Comfeel® Plus & Comfeel® Plus Transparent -  
Heals and Seals; Comforts and Protects

Protection and moist wound healing



## Key benefits

- Free from animal products, therefore suitable for those who don't use animal derivatives for religious or ethical reasons
- Smooth topfilm for low friction
- Soft and flexible for easy application and patient comfort
- Longer wear time<sup>1</sup> which reduces dressing change for the patient, therefore cost effective
- Mapping grid for wound size monitoring

# Protection and monitoring of skin and wound

Comfeel dressings are indicated for leg ulcers, pressure ulcers, post-operative wounds, superficial burns, superficial partial thickness burns, abrasions, donor sites and skin tears

## Comfeel® Plus Transparent

Protection and visual monitoring of skin and low exuding, non-infected wounds

- Transparency for easy wound monitoring
- Flexible and thin for easy application in difficult to dress areas
- Soft and discreet for patient comfort



Easy wound monitoring without wound disturbance

## Comfeel® Plus

Protection and treatment of low-to-moderate exuding, non-infected wounds

- Contains calcium alginate for enhanced absorption capacity<sup>2</sup>
- Bevelled edges for close adhesion to reduce risk of rolling edges<sup>3</sup>
- Soft and discreet for patient comfort



Wound size mapping and dressing change indication (gelling) help decide time for change

# Ensured protection. Faster healing<sup>4</sup>. Heals and Seals; Comforts and Protects

The Comfeel range of sizes and shapes delivers treatment flexibility for protection and monitoring of skin and wound

### Comfeel Plus Transparent



Size in cm	Item no.	Units per box	NHS	PIP
5x7	3530	10	ELM036	236-7332
10x10	3533	10	ELM060	282-7780
15x15	3539	5	ELM037	292-7798
20x20	3545	5	ELM190	282-7806

### Comfeel Plus



Size in cm	Item no.	Units per box	NHS	PIP
4x6	3146	30	ELM064	285-2945
10x10	3110	10	ELM351	222-9904
15x15	3115	5	ELM353	222-9912
20x20	3120	5	ELM164	222-9920

### Comfeel Plus Transparent



Size in cm	Item no.	Units per box	NHS	PIP
5x15	3547	10	ELM096	301-0618
5x25	3548	5	ELM163	301-0600
9x14	3536	10	ELM002	236-7340
9x25	3537	5	ELM189	296-2561
15x20	3542	5	ELM004	236-7324

### Comfeel Plus Contour



Size in cm	Item no.	Units per box	NHS	PIP
6x8	3280	10	ELM033	228-4818
9x11	3283	10	ELM034	228-4826

### Comfeel Plus Sacral



Size in cm	Item no.	Units per box	NHS	PIP
18x20	3285	5	ELM001	252-9055

The Comfeel range does not contain pectin, known allergens such as latex or colophony, or animal derivatives such as gelatin

1. Holm C, Petersen JS, Grønbæk F et al. Effects of occlusive and conventional gauze dressings on incisional healing after abdominal operations. Eur J Surg 1998. 164:179-183  
 2. Thomas S et al. An in-vitro comparison of the physical characteristics of hydrocolloids, hydrogels, foams, and alginate/cmc fibrous dressings. SMTL 2004  
 3. Ann Goodhead, Clinical efficacy of Comfeel Plus Transparent Dressing, British Journal of Nursing, 2002, Vol 11, No4, pp284-287)  
 4. Singh A, Halder S, Menon GR et al. Meta-analysis of randomized controlled trials on hydrocolloid occlusive dressing versus conventional gauze dressing in the healing of chronic wounds. Asian J Surgery 2003. 27 (4):326-332  
 \* Bacteria and virus impermeability: Coloplast data on file

