

Physiotulle®

- A non-stick hydrocolloid based wound contact layer



Physiotulle®


A non-stick hydrocolloid based wound contact layer

Comfort and Protection:

- **Pain and trauma free removal:** The dressing prevents adherence to the wound bed to ensure patient comfort, and can be removed without damaging newly formed granulation tissue
- **Thin and flexible:** Conforms to wound bed giving intimate contact
- **Easy to handle and use:** Does not stick to gloves
- **Hydrocolloid based contact layer:** Helps to support optimal wound conditions, so achieving faster wound healing¹
- **Wear time:** 2 days depending on wound type, exudate level & progression

What features of Physiotulle will benefit my patients?

- Upon contact with wound exudate, the hydrocolloid particles and the petrolatum maintain a moist environment beneficial for wound healing
- Physiotulle's thin, flexible structure conforms to the wound bed and is easily removed
- The petrolatum and hydrocolloid particles ensure that the dressing does not dry out or stick to the wound bed, so enhancing patient comfort



On what types of wounds would I use Physiotulle?

Physiotulle is indicated for the treatment of acute and chronic wounds such as:

- Superficial partial thickness burns
- Donor sites
- Post-operative wounds
- Skin abrasions
- Pressure ulcers
- Leg ulcers
- Patients with dry or fragile skin



Physiotulle Non-Adherent Wound Contact Layer

	Size (cm)	Qty	Code	NHS	PIP
	10x10	10	3910	EKB061	292-7986
	15x20	10	3915	-	292-8034

Ostomy Care
Urology & Continence Care
Wound & Skin Care

Coloplast develops products and services that make life easier for people with very personal and private medical conditions. Working closely with the people who use our products, we create solutions that are sensitive to their special needs. We call this intimate healthcare. Our business includes ostomy care, urology and continence care and wound and skin care. We operate globally and employ more than 7,000 people.

The Coloplast logo is a registered trademark of Coloplast A/S. © 09.2011
All rights reserved Coloplast A/S, 3050 Humlebæk, Denmark.

References
1. Humbert, et al. Acute and Superficial Chronic Wound Management in Outpatient Care. Impact Médecine, France – March 2003.



5 024683 009827 >



Coloplast Ltd
Peterborough Business Park
Lynchwood
Peterborough
PE2 6FX

www.coloplast.co.uk/woundandskincare

CM620N